
THUASNE

Ligaflex® Immo 0°

- fr** Attelle d'immobilisation du genou à 0°
en Knee immobilisation splint at 0°
de Schiene zur Ruhigstellung des Knies in einem 0°-Winkel
nl Immobilisatiebrace voor de knie op 0°
it Tutore di immobilizzazione per il ginocchio a 0°
es Férula de inmovilización de la rodilla a 0°
pt Tala para imobilização do joelho em 0°
da Immobiliserende knæskinne ved 0°
fi Immobilisoiva polvituki, taivutus 0°
sv Immobiliseringsskena för knät 0°
el Νάρθηκας ακινητοποίησης γονάτου σε κάμψη 0°
cs Imobilizační ortéza kolene 0°
pl Orzeza do unieruchomienia kolana pod kątem 0°
lv Celgala imobilizācijas sāite 0°
lt Kelio imobilizavimo 0° padėtyje įtvartas
et Põlve fikseerimise lahas 0° juures
sl Opornica za imobilizacijo kolena s kotom 0°
sk Imobilizačná kolenná ortéza pri 0° uhle
hu Oldalsínes térdrögzítő 0°-os szög
bg Шина за имобилизация на коляно на 0°
ro Atelă pentru imobilizarea genunchiului la 0°
ru Шина для имобилизации коленного сустава под углом 0°
hr Udlaga za imobilizaciju koljena, 0°
zh 0° 膝关节固定夹板
ar جبيرة لتثبيت الركبة على زاوية صفر درجة




		< 150 cm	150 - 175 cm	> 175 cm
	S	35 - 45 cm	 40 cm	x
	M	45 - 52 cm	x	 50 cm
	L	53 - 62 cm	x	 60 cm
	XL	63 - 75 cm	x	 60 cm

THUASNE
120, rue Marius AUFAN
92300 Levallois-Perret - France

		< 150 cm	150 - 175 cm	> 175 cm	
	S	35 - 45 cm		40 cm	
	M	45 - 52 cm	x		
	L	53 - 62 cm			
	XL	63 - 75 cm			
				50 cm	
					60 cm

fr	Taille patient	Circonférence cuisse	Hauteur
en	Patient height	Thigh circumference	Height
de	Größe des Patienten	Oberschenkelumfang	Höhe
nl	Maat van de patiënt	Dijomtrek	Hoogte
it	Altezza del paziente	Circonferenza coscia	Lunghezza
es	Tamaño del paciente	Circunferencia del muslo	Altura
pt	Tamanho do paciente	Circunferência da coxa	Altura
da	Patientens højde	Lårets omkreds	Højde
fi	Potilaan koko	Reiden ympärysmitta	Korkeus
sv	Patientens storlek	Omkrets runt lår	Höjd
el	Μέγεθος ασθενούς	Περιφέρεια μηρού	Ύψος
cs	Výška pacienta	Obvod stehna	Délka
pl	Rozmiar pacjenta	Obwód uda	Wysokość
lv	Pacienta augums	Augstīlba apkārtmērs	Augstums
lt	Paciento dydis	Šlaunies apimtis	Aukštis
et	Patsiendi suurus	Reie ümbermõõt	Kõrgus
sl	Velikost pacienta	Obseg stegna	Višina
sk	Veľkosť pacienta	Obvod stehna	Výška
hu	A beteg mérete	Comb körmérete	Magasság
bg	Размер на пациента	Бедрена обиколка	Височина
ro	Talia pacientului	Circumferința coapsei	Înălțime
ru	Размер пациента	Окружность бедра	Высота
hr	Veličina pacijenta	Promjer bedra	Visina
zh	患者身高	大腿圍	高度
ar	مقاس المريض	محيط الفخذ	الارتفاع

Immobilisation
Immobilisation
Ruhigstellung
Immobilisatie
Immobilizzazione
Inmovilización
Imobilização
Støtte
Immobilisoi
Immobilisering
Ακινητοποίηση
Znehybnění
Unieruchomienie
Imobilizācija
Imobilizavimas
Liikumatuks muutmīne
Imobilizacija
Imobilizácia
Rögzítés
Обездвижане
Imobilizare
Имобилизация
Imobilizacija
固定
التثبيت

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Conserver cette notice.

en

KNEE IMMOBILISATION SPLINTS AT 0° AND 20°

Description/Destination

The device is intended only for the treatment of the indications listed and for patients whose measurements correspond to the sizing table.

The device is composed of :

- 3 removable, radiolucent, rigid stays,
- 4 self-fastening closure and support straps,
- 1 removable anti-slip device,
- Foam padding for patient comfort.

Bilateral product.

Available in 1 size for the 40 cm high version (0° version only), 4 sizes for the 50 cm high version and 3 sizes for the 60 cm high version.

Composition

Textile components: Polyester-Polyurethane-Cotton-Polypropylene-Polyamide- Synthetic styrene-butadiene rubber (SBR).
Rigid components: Aluminium - Polyamide.

Properties/Mode of action

Knee joint stabilisation via the rigid stays.

Total immobilisation at 0° or 20°.

Adjustable closure and support via self-fastening straps.

Indications

Post-operative immobilisation.

Post-traumatic immobilisation.

Knee pain management.

Contraindications

Do not use the product if the diagnosis has not been confirmed. Do not use in the event of known allergy to any of the components. Do not apply the product in direct contact with broken skin or an open wound without an adequate dressing.

History of venous or lymphatic disorders.

Do not use in the event of major venous thromboembolic history without thromboprolyaxis.

Do not use in the event of unstable fractures, fractures of the lower part of the tibia or fibula, or fractures of the upper part of the femur.

Precautions

Verify the product’s integrity before every use.

Do not use the device if it is damaged.

Choose the appropriate size to fit the patient, referring to the size chart.

It is recommended that a health care professional supervises the first application.

The use of a medical device by a child, should be done under the supervision of an adult or a healthcare professional.

Strictly comply with your healthcare professional’s prescription and recommendations for use.

This product is intended for the treatment of a given condition. Its duration of use is limited to this treatment only.

It is recommended that the venous thromboembolic risk be assessed before any immobilisation is indicated. Follow the advice of a healthcare professional.

For hygiene and performance reasons, do not re-use the product for another patient.

It is recommended to adequately tighten the device to achieve support/immobilisation without restricting blood circulation.

In the event of discomfort, significant hindrance, pain, variation in limb volume, abnormal sensations or change in colour of the extremities, remove the device and consult a health care professional.

In the event of any modification in the product’s performance, remove it and consult a healthcare professional.

Prolonged immobilisation can sometimes cause muscle weakness.

Do not wear the product when driving a vehicle.

Do not use the device in case of application of certain products on the skin (creams, ointments, oils, gels, patches…).

Undesirable side-effect

This device can cause skin reactions (redness, itching, burns, blisters…) or wounds of various degrees of severity.

Any serious incidents occurring related to the device should be reported to the manufacturer and to the competent authority of the Member State in which the user and/or patient is resident.

Instructions for use/Application

Choose the splint height, then the size.

Measurements should be taken at the widest part of the thigh.

It is recommended that a health care professional supervises the first application.

It is recommended to wear the device directly on the skin, unless contraindicated.

- Open the splint and lay it flat open ❶.
- Place the opened splint under the leg ❷, the wider edge towards the top of the leg.

The patella must be centred in the cut-away area.

N.B.: The top of the splint must not impede flexion of the leg at the groin.

- Shaping the stays: if needed, it is possible to shape the stays in order to avoid potential pressure points or the splint slipping. Apply a slight pressure in the desired direction.

Repeat the process on the other stays if needed.

If subsequent shaping of the stays is needed, consult a healthcare professional.

- Engage the self-fastening straps in the buckles. Tighten the straps in the order described (A,B,C,D) ❸. Once all the straps are fastened, readjust the tightening.

If needed, reduce the strap length by cutting the edges with scissors.

Care/Maintenance

Product can be washed in accordance with the conditions shown on this instruction leaflet and on the label. Hand wash. Close the self-fastening tabs before washing. Remove stays before washing. Stays can be washed with a moist cloth. Do not use detergents, fabric softeners or aggressive products (products containing chlorine). Squeeze out excess water. Dry away from any direct heat source (radiator, sun, etc.). Dry flat.

Storage

Store at room temperature, preferably in the original packaging.

Disposal

Dispose of in accordance with local regulations.

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