

THUASNE

LombaStab vario

- (fr) Ceinture de soutien lombaire avec coussin lombaire - "Serrage facile et précis"
- (en) Lumbar support belt with lumbar cushion - "Easy, precise tightening"
- (de) LWS-Stabilisierungsohrthese mit Pelotte und Zugsystem - „einfaches und präzises Justieren“
- (nl) Lumbale rugbrace met lendensteun - "Eenvoudige en nauwkeurige aanspanning"
- (it) Cintura di sostegno lombare con cuscino lombare - "Serraggio facile e preciso"
- (es) Faja de sujeción lumbar con almohadilla lumbar de fijación fácil y precisa
- (pt) Cinta de suporte lombar com almofada lombar - "Aperto fácil e preciso"
- (da) Rygstøttebælte med lændepude - "Nem og præcis stramning"
- (fi) Selkä- ja lantiotukivyyö lannetyynyillä - Helposti ja tarkasti kiristettävä
- (sv) Ryggstöd för ryggen med kudde för ländryggen - "enkel och precis åtdragning"
- (no) Ortose med lendepute for stabilisering av korsryggen - "Enkel og presis stramming"
- (el) Ζώνη οσφυϊκής υποστήριξης με οσφυϊκό μαξιλάρι - "Εύκολη και ακριβής σύσφιξη"
- (cs) Bederní podpurný pás s bederním polštářkem - „Jednoduché a přesné dotáhnutí“
- (pl) Ortezta łądźwiowo-krzyżowa z poduszką łądźwiową - „Łatwe i precyzyjne mocowanie“
- (lv) Muguras jostas daļas atbalsta josta ar polsterējumu - „Viegla un precīza savilkšana“
- (lt) Palaikomasis juosmens diržas su juosmens pagalvėle - „Lengvai ir tiksliai suveržiamas“
- (et) Selja-ja nimmepiirkonna tugivöö nimmepadjaga - „Lihtne ja täpne pingutus“
- (sl) Oporni ledveni pas z ledveno blazino - „Enostavno in natančno vpetje“
- (sk) Podporný driekový pás s odnímateľnou driekovou pelotou - „Jednoduché a rovnomerné utiahnutie“
- (hu) Ágyéki fűző deréktámasztó párnával - „Egyszerű és pontos rögzítés“
- (bg) Лумбален колан с лумбална възглавница - „Лесно и точно затягане“
- (ro) Centură de susținere lombară cu pernă lombară - „Strângere ușoară și precisă“
- (ru) Опорный поясничный корсет с поясничной подушкой — «Легко и просто затянуть»
- (hr) Pojas s jastučićem za potporu lumbalnog dijela - „Jednostavno i precizno stezanje“
- (zh) 易于调节的带有腰部衬垫的腰部稳定护带
- (ja) 腰椎パッド付き腰サポートベルト - “簡単で正確な締め付け”
- (ar) حزام لدعم أسفل الظهر مزود بوسادة - "شد سهل ودقيق"

	3	4
	60 - 80 cm	1
	75 - 95 cm	2
	90 - 110 cm	3
	105 - 125 cm	4
	120 - 140 cm	5

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LUMBAR SUPPORT BELT WITH LUMBAR CUSHION – “EASY, PRECISE TIGHTENING”*

(en)

Description

- Elastic fabric ❶ (polyamide, polyester, elastane, viscose)
- Back panel lined with Coolmax fabric (polyamide, polyurethane)
- Four reshapable steel spring back stays (brushed tempered steel)
- Two flexible plastic front stays (polyoxymethylene)
- Pads and covers ❷ (polyamide)
- Brushed fleece panels (polyamide, polyurethane)
- Hooks (polyamide)
- Handles ❸ (polypropylene, SEBS)
- Laces (polyethylene)
- Lumbar cushion ❹ (polyamide)
- Cover (polyamide, cotton, polyester, elastane, acetate)

Properties / Mechanism of action

The LombaStab vario belt relieves low back pain, supporting the lumbar spine: compartmental lumbar unloading is therefore achieved via increased intra-abdominal pressure. This device uses back pads, implementing a pulley-like system to distribute pressure.

Indications

Non-specific low back pain (acute, subacute and chronic stages) in patients with a waist size of between 60 and 140 cm. Lumbar spinal stenosis.

Destination:

This brace should only be used for back treatment.

Instructions for use/Fitting instructions

- Choose the appropriate size to fit the patient, referring to the size chart (indicated on the box and on the instruction leaflet).
- Before use, ensure that the product is not damaged.
- Wear the belt on top of a t-shirt or thin vest.
- Open the belt and lay it out completely flat so that you can see the handles. Check that the two tightening handles ❷ are fully retracted. If necessary, unhook them and separate the two back parts by pulling the ends of the belt, then fasten the handles again.

Stabilisation phase

- If recommended by your healthcare professional, position the lumbar cushion ❹ in the belt using the touch fastening system (triangular design pointing down) ❶.
- Position the belt around the waist, placing your hands in the hand loops located on the front panels ❸. The tightening back pads ❷ must be positioned outwards (i.e. not against the skin). The lacing section must be centred over the spine: the pads should be positioned on either side of the spine. The bottom of the belt must stop at the top of the crease between the buttocks. The Thuasne logo indicates which way up the belt should be put on.
- Fasten the belt using the wide touch fastening straps: right side over left side, without pulling too tight or applying any particular pressure ❹.
- Adjust the tightness using the two handles: pull them forwards then fold them flat against the front of the belt, where they attach like a touch fastener ❸.
- These two handles can be used to adjust the level of compression as required throughout the day ❷.

Mobilisation phase

Remove lumbar cushion ❹ then use the belt in the same manner as during the immobilisation phase.

Contraindications

- Do not use on children and adolescents who have not finished growing, as no data are currently available.
- Do not apply the product in direct contact with skin (risk of pinching).
- Do not use on pregnant women.
- Do not use in the event of known allergy to any of the materials (see product composition in the Description section).
- Do not use in the event of a hiatus hernia.
- Do not use in the event of bone cancer with spinal metastases.
- Do not use in the event of deficient respiratory functions.
- Do not wear the product while sleeping.
- Do not wear the product in a medical imaging machine.

Precautions

It is recommended that a health professional adjust the stays to the shape of the patient's back before the belt is worn for the first time.

For temporary use, combined with functional rehabilitation if necessary, following the advice or prescription of a health professional.

In the event of discomfort or pain, consult your health professional.

Do not use the device if it is damaged.

Store the device at room temperature, in its original packaging. For hygiene and performance reasons, do not re-use the product for another patient.

For safety reasons, do not wear the belt close to machines (risk of becoming caught).

Dispose of the product with household waste (no special disposal conditions).

Unwanted side effects

Any serious incidents occurring related to the device should be reported to the manufacturer and to the competent authority of the Member State in which the user and/or patient is resident.

Washing instructions

- Fasten the touch fasteners before washing.
- Machine wash in a wash bag (belt and lumbar cushion cover) at a maximum temperature of 30 °C (thirty degrees Celsius).
- Clean the lumbar cushion with a damp cloth.
- Do not dry clean.
- Do not wash the device using aggressive products containing chlorine or oxygenated whitening products.
- Do not tumble dry.
- Do not iron
- Squeeze out excess water and dry flat away from heat sources.

*Study conducted in-house on a panel of 13 people, July 2019. This medical device is a regulated health product with CE marking.

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Thuasne Benelux B.V.
Thuasne Scandinavia AB
Thuasne UK Ltd
Thuasne (Shanghai) Medical Devices Co., Ltd.