THUASNE

LigaStep Pedi Walker

fr	Botte de marche pédiatrique à montants latéraux rigides
en	Pediatric walking boot with rigid side uprights
de	Pädiatrische Unterschenkel-Fußorthese mit steifen Seitenstützen
nl	Walker voor kinderen met laterale, rigide staanders
it	Stivale da deambulazione pediatrico con scocca laterale rigida
es	Bota de caminar pediátrica con montantes laterales rígidos
pt	Bota de caminhada pediátrica com suportes laterais rígidos
da	Pædiatrisk støttestøvle med stive stivere i siden
fi	Kokonaan jääykkärunkoinen lasten kävelysaapas
sv	Pediatriskt fotstöd med styvt lateralt skal
el	Παιδιατρικός ύάρθηκας ποδοκνημικής με άκαμπτες πλευρικές ράβδους
cs	Dětská fixační ortéza s pevnými bočními výztuhami
pl	But ortopedyczny pediatryczny ze sztywnymi wzmocnieniami boczonymi
lv	Stingrs sānu balsta bērnu pastaigu zābaks
lt	Eiti padedantis įtvaras vaikams su standžiais šoniniais statramsčiais
et	Jäikade külgededega pediaatriline tugisaabas
sl	Pediatrični škorenj s trdima stranskima opornikoma
sk	Detská členková dlaha s pevnými postrannými dlahami
hu	Gyerek járótalpas bokarögzítő oldalmerevitővel
bg	Педиатрична ортеза за ходене mun „бомуш“ с твърди страннични опори
ro	Cizmă pediatrică cu montanți lateralii rigizi
hr	Pedijatrijska čizma sa čvrstim kalupom
ru	Детский сапог для ходьбы с жесткими боковыми стойками
zh	儿童助行靴，带坚固护腿束带
ar	حذاء للمشي للأطفال ذو دعامتين جانبية صلبة



S	21 - 25
M	25 - 31
L	31 - 34

THUASNE
120, rue Marius Aufan
92300 Levallois-Perret - France

www.thuasne.com

© Thuasne - 6162501 (2020-10). ©Studio Caterin

3	4
S	21 - 25
M	25 - 31
L	31 - 34

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Conserver cette notice.

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Liner fully removable for washing.
Hand wash the liner in cold water.
Wash the rigid part with a moist cloth.
Squeeze out excess water. Dry away from any heat source (radiator, sun, etc.).

Do not use detergents, fabric softeners or aggressive products (products containing chlorine).

Storage

Store at room temperature, preferably in the original packaging.

Disposal

Dispose of in accordance with local regulations.

Keep this instruction leaflet.

PEDIATRIC WALKING BOOT WITH RIGID SIDE UPRIGHTS

Description/Destination

Walking boot with rigid side uprights to immobilize and stabilize the lower third of the leg, ankle and foot.

The device is intended only for the treatment of the indications listed and for patients whose measurements correspond to the sizing table.

Composition

Rigid components: polyamide, fiberglass, rubber, ethylene-vinyl acetate.

Fabric components: polyamide, polyurethane, polyester.

Properties/Mode of action

Immobilizes and stabilizes the lower third of the leg, ankle and foot.

90° rigid uprights on both sides

Designed to ensure minimal dorsiflexion of the metatarsal bone.

Curved, anti-slip sole that offers patients a natural, comfortable walking motion.

Product can be worn on either foot.

Indications

Ankle/foot post-operative/post-traumatic immobilization.

Severe ankle/foot sprain.

Stable fracture of the lower third of the leg (tall version).

Conservative treatment of ankle/foot fracture.

Contraindications

Do not use the product if the diagnosis has not been confirmed.

Do not apply the product in direct contact with broken skin.

Do not use in the event of known allergy to any of the components.

Do not use in the event of unstable fractures or fractures of the upper part of the tibia or fibula.

Do not use for patients weighing > 68 kg.

History of venous or lymphatic disorders.

Precautions

Verify the product's integrity before every use. Do not use the device if it is damaged.

Choose the appropriate size to fit the patient, referring to the size chart.

It is recommended that a healthcare professional supervises the first application.

Strictly comply with your healthcare professional's prescription and recommendations for use.

The systematic use of a long sock is recommended when wearing the device.

Ensure correct position and avoid folds when placing foam liner in the device.

Check the condition of the affected limb daily (with particular attention for patients with sensory deficit). Do not bring the product into direct contact with greasy substances (ointment, cream, etc.).

Do not bring the product into direct contact with a dressing.

In the event of discomfort, significant hindrance, pain, variation in limb volume, abnormal sensations or change in colour of the extremities, remove the device and consult a healthcare professional.

Do not wear the product in a medical imaging device. Do not expose the product to extreme temperatures (> 60°C).

It is recommended to adequately tighten the device to achieve immobilisation without restricting blood circulation. In the event of pain, tingling or numbness, loosen the straps or deflate the air pads and ask a healthcare professional for advice.

For hygiene and performance reasons, do not re-use the product for another patient.

If the device comes into contact with water, dry the textile part and wipe the rigid part well with a dry cloth.

Prolonged immobilisation can sometimes cause muscle weakness.

When resuming walking, be careful of risks of falling.

Garment care

Product can be washed in accordance with the conditions shown on this instruction leaflet and on the label.

Undesirable side-effect

This device can cause skin reactions (redness, itching, burns, blisters...) or wounds of various degrees of severity.

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