
THUASNE®


SLEEQ® MAX TLSO

en	Postoperative/post-injury TLSO (thoracolumbar sacral orthosis).....	6
de	TLSO (Thorakal-Lumbal-Sakral-Orthese) zur postoperativen/posttraumatischen Versorgung.....	7
fr	TLSO (orthèse thoraco-lombo-sacrée)postopératoire/post-traumatique	8
nl	TLSO (thoraco-lumbo-sacrale orthese) postoperatief/posttraumatisch	9
it	TLSO (ortesi toraco-lombo-sacrale) postoperatoria/post-traumatica	10
es	TLSO (órtesis toracolumbosacra) posoperatoria/postrumática.....	12
pt	TLSO (ortótese toracolumbosacral)pós-operatório/pós-traumático	13
da	Postoperativ/posttraumatisk TLSO (torakolumbalsakral ortose).....	14
fi	TLSO (rintakehän, lannerangan ja ristiselän tuki) leikkauksen/trauman jälkeen.....	15
sv	TLSO för post-op/post-skada = Ortos för bröst-, länd- och sakralområdet.....	16
el	Θωρακο-οσφυο-ιερός κηδεμόνας (TLSO) μετεγχειρητικός / μετατραυματικός.....	18
cs	TLSO (torakolumbosakrální ortéza) pooperační/posttraumatická.....	19
pl	TLSO (orteza piersiowo-lędźwiowo-krzyżowa) pooperycyjna/pourazowa	20
lv	TLSO (krūšu-jostas-krustu daļas ortoze) lietošanai pēc operācijas / pēc traumas	21
lt	Po operacijos ir (arba) po traumos dėvimas TLSO (krūtinės-juosmens-kryžkaulio įtvartas).....	23
et	TLSO (rindkere- nimm- ristluu ortoos) operatsiooni- /traumajärgseks kasutuseks.....	24
sl	TLSO (prsna-ledvena-krizna opornica) postoperativna/posttravmatska.....	25
sk	TLSO (torakolumbosakrálna ortéza) pooperačná/posttraumatická.....	26
hu	TLSO (mellkasi-ágyéki-keresztscsonti ortézis) műtét utáni /poszttraumás.....	27
bg	TLSO (торакo-лумбо-сакрална ортеза) – постоперативна /посттравматична	28
ro	OTLS (orteză toraco-lombosacrală) postoperatorie/posttraumatică.....	30
ru	TLSO = торакoлyмбoсакральный ортез послеоперационный/посттравматический	31
hr	TLSO (potporno rasteretna ortoza za kralješnicu) postoperativna/posttraumatska	32
zh	术后/创伤后 TLSO (胸腰骶矫形器)	33
ar	TLSO (دعامة صدرية قطنية عجزية) لفترة ما بعد الجراحة/ما بعد الإصابة بالرضوح.....	35


1.0


2.0


3.0


4.1

4.2


4.3


* encoche sternale - Brustbeineinkerbung - sternale inkeping - appoggio sternale - escotadura esternal - entalhe do esterno - sternalt hak - rintalastan lovi - bröstskorg - στερνική εντομή - hrudní zářez - wcięcie mostka - krūšu kaula bedrīte - krūtinkaulio viršus - rinnakuosa sälguga - suprasternalna ząreza - hrudný zářez - szegycsoni bevágás - зр'ядна вгальбнатина - fosă sternală - грудная выемка - prsna kost - 胸骨切迹 - التلمة القصية

5.0

5.1


5.3


* ou - oder - of - o - o - ou - eller - tai - eller - n - nebo - lub - vai - arba - vði - ali - alebo - vagy - или - sau - или - ili - 或 - 或

5.2


5.4


6.0

6.1

6.2

en	Universal size
de	Einheitsgröße
fr	Taille unique
nl	Eén maat
it	Taglia unica
es	Talla única
pt	Tamanho único
da	Kun én størrelse
fi	Yksi koko
sv	En storlek
el	Ένα μέγεθος
cs	Univerzální velikost
pl	Rozmiar uniwersalny
lv	Universāls izmērs
lt	Vienas dydis
et	Üks suurus
sl	Ena velikost
sk	Univerzálna veľkosť
hu	Egy méret
bg	Един размер
ro	Mărime universală
ru	Единый размер
hr	Univerzalna veličina
zh	均碼
ar	مقاس وحيد

Immobilisation
Ruhigstellung
Immobilisation
Immobilisatie
Immobilizzazione
Inmovilización
Imobilização
Støtte
Immobilisoi
Immobilisering
Ακινητοποίηση
Imobilizace
Unieruchomienie
Imobilizācija
Imobilizavimas
immobiliseeriv
Imobilizacija
Imobilizácia
Rögzítés
Обездвижване
Imobilizare
Иммобилизация
Imobilizacija
固定
التثبيت

FDA
REGISTERED

POSTOPERATIVE/POST-INJURY TLSO (THORACOLUMBAR SACRAL ORTHOSIS)

Description/Use

Spinal orthosis for active relief and correction of the lumbar and the thoracic spine in the sagittal plane.

Rigid tri-planar control TLSO. Universal size and customizable features to provide support for the patient.

Composition

Textile components: polyamide - polyester - synthetic rubber (latex-free).

Rigid components: polypropylene - stainless steel - brass - aluminum.

Properties/Mode of action

Rigid tri-planar immobilisation shell design featuring customizable structural components and an adjustable sternal bar.

Adjustable compression is delivered via pull-tabs connected to a corset-style lacing system.

Indications

Sagittal/coronal/transverse plane motion restriction.

Stable vertebral fracture.

Increased intracavitary pressure to reduce spine and disc loading.

Mechanical or discogenic thoracic/lumbar and extremity pain.

Postoperative spinal stabilisation/fusion.

Degenerative disc disease.

Spondylolysis.

Post-injury to facilitate healing.

Spinal stenosis.

Spondylolistesis.

Spondylolysis.

Laminectomy.

Discectomy.

Facet syndrome.

Spinal muscle weakness.

Arthritis.

Low back sprains and strains.

Prophylactic back support.

Contraindications

Do not use the product if the diagnosis has not been confirmed.

Do not use in the event of known allergy to any of the components.

Do not apply the product in direct contact with broken skin or an open wound without an adequate dressing.

Do not use in the event of circulatory, pulmonary or cardiovascular problems in patients for whom an increase in arterial blood pressure is not recommended.

Do not use on pregnant women.

Precautions

Additional accessories/spare parts are available to order.

Do not use the device in case of application of certain products on the skin (creams, ointments, oils, gels, patches...).

It is recommended that your healthcare professional supervises the first application.

Verify the product's integrity before every use.

Wear the device above thin clothing.

Do not use the device if it is damaged.

Choose the appropriate size to fit the patient, referring to the size chart.

Strictly comply with your healthcare professional's prescription and recommendations for use.

For hygiene and performance reasons, do not re-use the product for another patient.

Do not wear the product in direct contact with skin (risk of pinching).

Do not wear the device close to machines (risk of becoming caught).

Do not wear the device in a medical imaging machine.

Do not wear the product when driving a vehicle.

Do not wear the product while sleeping.

It is recommended to adequately tighten the device to achieve support without excessive compression.

In the event of discomfort, significant hindrance, pain, abnormal sensations, remove the device and consult a healthcare professional.

Undesirable side-effect

Any serious incident arising that is linked to the product should form the subject of a report to the manufacturers and the competent authority for the Member State in which the user and/or patient is domiciled.

This device can cause skin reactions (redness, itching, irritations...).

Instructions for use/Application

Preparation of the device by a healthcare professional:

1.0 Brace components

A Anterior panel assembly

B Lateral panels (x2)

C Wing shells (x2)

D Pull-tabs

E Crescent connectors

F Strap clip

G Sizing arcs

H Extension panel (*sold separately)

I Sizing guide

J Wing

K Sternal pad

L Sternal bar clip

M Sternal bar

N Strap anchor

2.0 Sizing

2.1 Prepare to fit:

Press the button and remove the straps from the strap clips on the sternal pad **K** **I**.

Remove the anterior panel **A**.

2.2 Adjust circumference:

Measure waist and see Sizing guide **I** **W**.

3.0 Attaching wing shells and lateral panels

3.1 Align wing shell **C** sew line arcs to sew line arcs on the wing.

Note: for size 4 use posterior arc near crescent connectors **E** on the back panel.

3.2 Tuck wing shell **C** inside the back panel.

3.3 Secure tab (A) to wing inside the back panel.

Repeat on opposite side.

3.4 Align the lateral panel **B**, then attach by folding over both tabs **B**.

Attach by folding over both tabs (B).

4.0 Trying on and fitting

4.1 Position the brace at coccyx level (I).

Note: have the patient hold the wings.

4.2 Remove the anterior panel **A**.

Wrap wings (I) and (2) attach to anterior panel.

4.3 Tighten the brace:

Move the pull-tabs **D** anteriorly and re-attach them on the wing.

4.4 Adjust the sternal bar:

Open the sternal bar clip **L** to adjust height, then close the clip to lock.

Sternal pad **M** should be located between 1 and 2 inches, or 2.5 to 5.0 cm (2 finger widths), below the sternal notch.

Sternal bar **M** can be bent to achieve proper fit with bending irons or over the edge of a table.

Note: only bend sternal bar above toothed area.

5.0 Modifications

5.1 The back panel, anterior panel **A** (shown), lateral panels and wing shells can be trimmed with scissors **(1)** and heat modified **(2)** to adjust lordosis, add scapula reliefs or create flares.

The steel sternal bar **M** can be bent to match anatomical contours and the wings can be attached to the crescent connectors **E** at up to 5 degree angles up or down to accommodate various waist to hip ratios.

Do not trim top center tabs **(3)** on the anterior panel.

5.2 Back panel removal for modification:

1 Open cover at the base of the back panel.

2 Remove the lumbar assembly.

3 Remove the plastic panel.

5.3 Adjust for over the shoulder or under axilla strapping option:

Detach the strap anchors **N** and reposition up for over the shoulder (A) or down for under axilla strap placement (B).

5.4 Removing or replacing the sternal bar (with available long bar option):

Remove anterior panel **A** from the wing.

Release sternal bar clip **L** and lower the bar.

Remove retainer clip (A) by pushing through bar while bending bar away from panel.

Pull sternal bar **M** up to remove from bar slot.

Insert new short sternal bar, (or optional long sternal bar) into bar slot and replace retainer clip (A).

Note: clip (A) must always be placed on the panel side of bar.

6.0 Removal

6.1 When removing the brace, always replace pull-tabs **D** to starting position.

6.2 Lay brace on a flat surface then pull both wings **1** outward to reset the tightening mechanism.

Care

The product can be washed in accordance with the conditions shown on this instruction leaflet and on the label. See product label. Hand wash. Wipe with a damp cloth. Do not tumble-dry. Do not iron. Squeeze out any excess water. Dry away from any direct heat source (radiator, sun, etc.). Do not expose the product to extreme temperatures. If the device is exposed to seawater or chlorinated water, make sure to rinse it in clear water and dry it.

Storage

Store at room temperature, preferably in the original packaging.

Disposal

Dispose of in accordance with local regulations.

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©Thuasne - CO0060000001 (2026-04)
THUASNE UK Ltd
Thuasne SAS - SIREN/RCS Nanterre 542 091 186